wwellion®

INOVO

FIA Analyzer

HANDBUCH
USER MANUAL
УКАЗАНИЕ ЗА УПОТРЕБА
MANUÁL
MANUEL D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
KORISNIČKI PRIRUČNIK
MANUALE
MANUAL DE INSTRUÇÕES
ANVÄNDARMANUAL
UPORABNIŠKI PRIROČNIK
PRÍRUČKA PRE POUŽÍVATEĽA

HANDBUCH

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USER MANUAL

Dear Wellion® INOVO FIA Analyzer User,

Thank you for choosing the Wellion® INOVO FIA Analyzer! This analyzer has the advantages of easy operation, strong stability, reliable results and simple maintenance. It can be used widely in many areas, such as hospital, medical center, clinic, pharmacy, laboratory, etc.

This User Manual applies to the Wellion® INOVO FIA Analyzer, it provides detailed descriptions of the analyzer's characteristics, graphical operation steps and other relevant information for your convenience. Please read this User Manual carefully before you use your analyzer and keep your User Manual in a safe place, you may want refer to it in the future.

Thank you again for choosing the Wellion® INOVO FIA Analyzer.

TABLE OF CONTENTS

Principle and Intended Use	34
Chapter 1: Understanding your analyzer	35
Basic Introduction	35
Important Safety Information	35
Warnings, Precautions and Limitations	35
Chapter 2: Components and structure	39
Components List	39
Structure	39
Chapter 3: Installation and Use conditions	41
Analyzer Installation	41
Use Conditions	41
Chapter 4: Operation instructions	43
Starting	43
Settings	43
UserMgr	47
Test	48
Batch Test	50
Result	51
Projects	53
Ratio	53
Chapter 5: Maintenance and Troubleshooting	55
Maintenance and Care	55
Troubleshooting	56
Chapter 6: Technical Information	58
System Specifications	58
Limitations	58
Index of symbols	59
Manufacturer's Warranty	60

Principle and Intended Use

The Wellion INOVO FIA Analyzer works with Wellion INOVO Test Kits as a system to measure multiple parameters like SARS-CoV-2 total antibody, SARS-CoV-2 Ag, etc. The analyzer is an immunofluorescence detection device that is based on photoelectric detection principle. With the fluorescence immunochromatography technology and tiny sample, the analyzer gives reliable result(s) in just few minutes, which is ideally suitable for professionals to make quick and correct decision at the point of care. It is intended for professional use outside the body (in vitro diagnostic use).

Basic Introduction

1. Product name: Wellion® INOVO FIA Analyzer
2. Basic size: 215 mm x 303 mm x 159 mm (W x D x H)
3. Total weight: 2.5 kg
4. Security classification of medical electrical equipment: type of protection against electric shock is Class I; pollution grade is Class II; facility category (overvoltage category) is Class II.

Important Safety Information

- For professional use only.
- For *in vitro* diagnostic use.
- Use only Wellion® INOVO Testkits with your analyzer. Do not use any other brand.
- Use the accessories provided by MED TRUST with your Wellion® INOVO FIA Analyzer. Do not use other accessories.
- Please read this User Manual carefully before using your analyzer and keep your User Manual in a safe place for referring to it in the future.
- Please comply with all the national standards for health and safety. This includes (but is not limited to) protective cover, gloves, and suitable experimental installation when operating and maintaining the analyzer.
- If the analyzer is used in a manner not specified by the manufacturer, the protection provided by the analyzer will be impaired.
- All products or objects, which come in contact with blood/urine/respiratory secretions, even after cleaning, should be handled as if capable of transmitting infectious diseases.

Warnings, Precautions and Limitations

Notices:

Warnings and precautions

(1) Warnings and precautions for operation place

Warnings:

- Please try not to use parallel socket so as to avoid overload and cause fire.
- A 24V/2.5A power adapter and an effectively grounded outlet must be used.
- Damaged or non-original or modified power cord has a risk of causing fire and electric shock. Do not over-bend or rolling the power cord so as to avoid a fire or electric shock.

- When the analyzer is loose or the parts are dropped out and damaged, please contact the manufacturer in time.
- Do not use this analyzer in unstable environments such as tilt, vibration, impact, etc.
- Do not place the device in a location where it is difficult to disconnect the device.
- Do not drop the analyzer or get it wet. No water or debris should enter into the Analyzer. If so, please contact the manufacturer.
- Do not take the analyzer apart. Taking the analyzer apart will void the warranty.
- Refer to the **Maintenance and Troubleshooting** for details on cleaning the analyzer.
- Keep the analyzer and accessories out of reach of children.

⚠ Notices:

- Turn off the power and unplug the power plug before moving the analyzer.
- Try to avoid vibration when moving the analyzer.
- Desktop supporting the analyzer should load greater than 2.5kg.
- The analyzer should be placed steadily. At least 5cm space should be saved to ensure air circulation and heat dissipation.
- The analyzer cannot be covered by any items.
- Do not place the analyzer near radio, television, printer, fax machines and other sources of interference.
- Do not use the analyzer with other instruments such as microwave and other high-frequency equipment.

(2) Warnings and precautions for using

⚠ Warnings:

- Read the User Manual carefully before starting the analyzer. The operator must receive professional training, and be familiar with the User Manual and operation method. The analyzer must be managed by a dedicated person.
- Set the test parameters under the guidance of the professionals.
- When handling potentially contagious substances such as animal specimens or reagents, protective gloves or other protective measures are required if exposed to skin.

⚠ Notices:

- Make sure the analyzer is in normal running status before use.
- Make sure that all wires are properly connected and secure.

- When using with other instruments together, be sure to read and clarify the operation precautions.
- After the test, confirm that the test cassette has been removed, the test cassette holder has been reset; and then turn off the power.

(3) Warnings and precautions for faults, storage and inspection

⚠ Warnings:

- Please stop the analyzer immediately when abnormal conditions occur. If the analyzer emits peculiar smell, smoke or liquid flows into the analyzer, please cut off the power immediately and contact your local dealer. Continuous use of the analyzer may cause fire, electric shock or personnel injury.
- In addition to the maintenance personnel from manufacturer and the ones authorized by manufacturer, others are not allowed to remove, modify or repair the analyzer. Any violation will lead to failure of normal warranty and maintenance. As the manufacturer, we will not bear any responsibility for possible personal injury and fire risk of electric shock caused by violation.

⚠ Notices:

- The analyzer and its parts need to be checked regularly. Please notify the manufacturer for repair or replacement if there are damage, cracks and other abnormal conditions.
- Use a clean soft cloth and non-corrosive cleanser to clean the analyzer surface to prevent scratching the shell and face plate of analyzer.

(4) Warnings and precautions for electromagnetic compatibility

⚠ Warnings:

- The analyzer is designed and tested for electric safety and compatibility. If used in home environment, it may cause radio interference and require taking protective measures.
- Do not use this analyzer near strong radiation sources such as unshielded RF sources. Otherwise it may interfere with proper operation of the analyzer.

⚠ Notices:

- The user should ensure that the analyzer is in an electromagnetic compatibility environment so that it can work properly.
- It is recommended to evaluate the electromagnetic environment before using the analyzer.
- This Analyzer complies with the noise immunity and emission requirements which is specified in the part of EN 61326-2-6.

(5) Limitation requirements for toxic and hazardous substances

- This analyzer meets the limitation requirements of toxic and hazardous substances according to Directive 2011/65/EU.

CHAPTER 2: COMPONENTS AND STRUCTURE

Please check whether any component is missing or damaged against the list below after opening the package.


Note: If some parts are missing or damaged, please contact your local dealer or manufacturer immediately.

Components List

Number	Component	Quantity
1	Wellion® INOVO FIA Analyzer	1
2	Power adapter	1
3	Power cable	1
4	Network cable	1
5	Printing paper	1
6	User Manual	1

Structure

Wellion® INOVO FIA Analyzer mainly consists of an optical unit, mechanical unit, control unit and display unit.

Analyzer (Figure 1)

(Figure 1)

Power adapter (Figure 2)

(Figure 2)

Printing Paper (Figure 3)

Recommended printing paper size :

Width: 57 mm


Diameter: no more than 29 mm

(Figure 3)

Code Chip (ID card) (Figure 4, additional purchase)


The code chip is a specific ID card with calibration curve, and can be only used with the test device. The code chip is supplied in the Wellion® INOVO Test Kit.

(Figure 4)

Test Device (Test cassette) (Figure 5, additional purchase)

The test device/Test cassette below is only supplied in the Wellion® INOVO Test Kit.

(Figure 5, additional purchase)

Analyzer Installation

Please use the analyzer under the proper conditions.

Unpacking and Checking

- Step 1:** Gently remove the analyzer and accessories from the packaging box. Save the packaging materials for future use.
- Step 2:** Check the analyzer and accessories if they are in good condition
- Step 3:** Connect the power adapter, and turn on the power switch to start the analyzer.

Notices:

- If there are any problems, please contact your local dealer or manufacturer immediately.

Installation and Debugging Procedures

- The analyzer should be placed in a clean and ventilated room with temperature between 10 °C and 30°C, relative humidity is less than 70%, away from direct sunlight.
- Make sure the vents are not obstructed, and that there is at least 5 cm of clearance around the analyzer.
- Connect the power adapter to the power interface of the analyzer and turn on the power.
- The place of the analyzer placed should be without strong vibration source and electromagnetic fields around.
- The analyzer has been debugged before shipping, and can be used directly.
- To ensure the analyzer works properly, please do not place any items on the analyzer at any time.

Use Conditions

Main Parameters

- Software version: Version 1
- Operation system: Linux
- Exciting light: LED or diode laser
- Connectivity: USB port, Ethernet port, double serial port (serial port 1: automatic LIS uploading; serial port 2: debugging)
- Screen: 7 inches 24 bits color touch screen
- Printer: Built-in thermal printer

Storage and Transportation Conditions

(1) Storage conditions

The analyzer should be stored in original packing box in a well-ventilated and clean room. The recommended environment temperature is - 20°C to 55°C. Relative humidity should be less than 93%. Harmful gas, flammable and explosive substances and corrosive gases are prohibited.

(2) Transportation conditions

- The analyzer with package is suitable for highway, railway, aviation and waterway transportation.
- Avoid severe vibration and shock during loading and transportation.
- Keep away from damp.
- Mix-package and mix-transportation with flammable and corrosive substances are prohibited.
- More specific requirements can be according to the contract.

Power Supply Voltage Requirements

Caution:


The power supply of the analyzer varies between 100-240V, 50/60Hz AC depending on the country where the analyzer is used. The input voltage is 24V DC. Please avoid short circuit and electric shock during use.

Starting

Turn on the analyzer, the analyzer will load and initialize the system (Figure 6, 7).


(Figure 6)

(Figure 7)

The following display will appear after initializing the system (Figure 8). Input the correct account name “admin” and password “admin”. Then click “Login” to log in the analyzer and to enter the main menu. Click “Keep” to keep the account name and password. Click “Auto Login” to log in the account automatically.

(Figure 8)


Settings

Please set up the analyzer according to your requirements before using. Click “Menu” button to enter the following interface (Figure 9).

(Figure 9)


Click the “Settings” icon to enter the following interface (Figure 10).

(Figure 10)

Institutions


Click the “Institutions” button to enter the following interface (Figure 11). You can check the institution name and institution address.

(Figure 11)

LIS/HIS


Click the “LIS/HIS” button to enter the following interface (Figure 12). To set the LIS uploading parameters, the upload-way needs to be selected at first. Then, the corresponding parameters be set.

(Figure 12)

Test settings

These settings include the instant and standard test settings. Click “Test settings” button to enter the following interface (Figure 13). You can set the serial number of test, length and alignment, etc. of test results` number.

(Figure 13)

1. “Serial:YMMDD-” to set the initial number of the test sample. „Serial: YMMDD -“ is the serial number formation rule. “YMMDD-“ represent the first six digits of the formed serial number according to the year, month and day of the test time.

Click the arrow of the “Serial subfix length” field to adjust the largest serial number of test; e.g. “3” will be shown as “001” and the maximum serial number will be “999”.

2. Input the number of days in the “Result expired time:” field to adjust the time to save the records. If you choose “0”, the analyzer will always save the record.

Slide the button “Serial alignment” to the right (green area) and turn on the function that the serial number of the test will be in a unified format, e.g. if length is set to “3” and function of “Serial alignment” is OFF, the sample`number will appear in the format “1”, if alignment is ON it will be “001”.

3. Slide the button “Sample No Auto Add” to the right to turn on the function that the sample ID will be added automatically when testing. Slide the button “Instant print” to the right to turn on the function of automatically printing the test result after Instant or STD test is finished. Slide the button “Auto test” to the right to turn on this function for the Instant Test. When the test device is inserted into the slot, the analyzer will move it in and measure the instant test cassette automatically and immediately.

Slide the button “Check card” to the right to turn on the function. The analyzer will check if the test cassette is inserted to the slot before a test is started.

4. RESULT TYPE „Both“ is used for the parameters need to display both

qualitative and quantitative results.

RESULT TYPE „Qualitative“ is used for the parameters need to display qualitative result only.

RESULT TYPE „Quantitative“ is used for the parameters need to display quantitative result only.

Note:


- Different parameter (/testkits) may require different result type. Make sure to choose the correct result type before test. Otherwise it may cause incorrect result.

Note:

- For all above listed functions and settings you need to click “Save” and “Confirm” to save the changes.
- For Batch test settings, refer to the Batch test interface.

System

Click “System” button to enter the following interface (Figure 14).


(Figure 14)

You can set the system date, time, turn on/off the touch sound, calibrate the screen, reset the analyzer and check the analyzer`slot function, in this interface.

1. Input the date and time, then click “Save” and “Confirm” to save the date and time.
2. Slide the “Touch beep” button to the right to turn on the sound for touching. Click “Save” and “Confirm” to save this setting.
3. Click “Calibrate” and “Confirm” to calibrate the touch screen.
4. Click “Reset” and “Confirm” to clear all data and restore factory settings.
5. Click “Check” to start the self-check of the test device slot.

About

Click “About” button to enter the following interface (Figure 15).


(Figure 15)

Article model, software release and full version information of the analyzer will be displayed.

Here, the user can also save the upgrade software package on an U disk, and then insert it into the USB port of the analyzer. The system will update after detecting the upgrading program.

UserMgr

Click the “UserMgr” button to enter the following interface (Figure 16).

(Figure 16)

1. Click “Logout” to log out the current account
2. Select Account and click “Modify” to adjust the user name and password (except for admin account). Click “Confirm” to save the setting.
3. Click “Add user” to add an account. Input the account, user name, password, and repeat the password. Choose the user type and click “Confirm” to add the account and sub-account.
4. Select account and click “Delete user” to delete the sub-account.

Test


Before you perform the test, you need to prepare the Wellion® INOVO FIA Analyzer and the compatible Wellion® Test Kit.

Operation Steps for Instant Test

Step 1: Connect the power adapter and turn on the power switch to start the analyzer.

Note: If you use the analyzer the first time or if the device has been moved to a different location/different environmental conditions, you have to wait 20 minutes, after turning on the analyzer, before making a measurement. 20 minutes is the preheating time for the analyzer to stable its system.

Step 2: Click “Test” button of main menu interface to enter the following interface (Figure 17).

(Figure 17)

Step 3: Remove the code chip (“ID card”) from the compatible Wellion® Test Kit, insert the code chip into the code chip (“ID card”) slot, and click the “Read ID” button. Then, confirm project/test item with the “Confirm” button. The analyzer will read the test item’s information automatically.

Note: If the analyzer displays error messages of project inconsistency („Cannot find matching Project“) after Reading the ID card, please check the Test settings interface, and turn off the „Scan barcode“ function in the Test settings.

Step 4: Select the “Sample type” and input the sample number “Sample No.,” if needed, after the code chip is recognized by the analyzer.

Step 5: Click “Info edit” button to input more detail information of the patient.

Step 6: Slide the Test type button to the left to select “Instant test”.

Step 7: Remove the test device (test cassette) from the Wellion® INOVO Test Kit and from the foil pouch, and refer to the test kit package insert for the specific test instructions (e.g. sample addition,

reaction/incubation time).

Step 8: After the incubation time is completed, insert the test device into the test device slot*, and then, click "Test". After pressing "Test", the analyzer will start the measurement immediately.

***Note:** The test cassette front has to be completely moved inside the analyzer slot, without gap at the front. Carefully check the right direction to move inside which is indicated by an arrow imprint on the test cassettes upper side.

Step 9: The analyzer will display the test result after measurement of the test cassette.

Insert this direction into the analyzer' cassette slot.

Operation Steps for Standard Test

Step 1 to 5 are the same as the operation steps for instant test.

Step 6: Slide the Test type button to the right to select "STD test" (=Standard Test).

Step 7: Remove the test device (test cassette) from the Wellion® INOVO Test Kit, and from the foil pouch, and refer to the test kit package insert for the specific test instructions (e.g. sample addition, reaction/incubation time).

Step 8: Insert the test device into the test device slot*, and then click "Test". The analyzer will start and display the incubation countdown now on the display and will incubate the test cassette inside the analyzer for the respective reaction time. Please wait for the necessary incubation time for each test before a measurement is done.

***Note:** The test cassette front has to be completely moved inside the analyzer slot, without gap at the front. Carefully check the right direction to move inside which is indicated by an arrow imprint on the test cassettes upper side.

Step 9: The analyzer will start to measure the test cassette automatically after the incubation time is completed, and, then, will display the test result.

Insert this direction into the analyzer' cassette slot.

Batch Test

Before you perform the test, you need to prepare the Wellion® INOVO FIA Analyzer and the Wellion® INOVO Testkit.


Step 1: Connect the power adapter and turn on the power switch to start the analyzer.

Note: If you use the analyzer the first time or if the device has been moved to a different location/different environmental conditions, you have to wait 20 minutes, after turning on the analyzer, before making a measurement. 20 minutes is the preheating time for the analyzer to stable its system.

Step 2: Enter the “Test” interface of the Main menu. Remove the code chip (“ID card”) from the compatible Wellion® Test Kit, insert the code chip into the code chip (ID card) slot, and click the “Read ID” button. Then, confirm project/test item with the “Confirm” button. The analyzer will read the test item’s information automatically. Click „Menu“ to return to the main menu.

Note: If the analyzer displays error messages of project inconsistency (“Cannot find matching Project”) after Reading the ID card, please check the Test settings interface, and turn off the “Scan barcode” function in the Test settings.

Step 3: Click the “Batch Test” icon of the main menu interface to enter the following interface (Figure 18).

(Figure 18)

Step 4: Select the “Sample type” and, then, check/select the test item “Project” to define the incubation time. The “Timer” in the screen will display the incubation time.

Step 5: Click the “Add” button to define and to increase the number of samples you will test. The analyzer will generate the sample number “Sample No.” automatically after the sample is added to this platform.

Note: If user wants to delete the redundant sample, select the respective sample and click the “Delete” and “Confirm” button. To modify the

Sample No., select the respective sample and click “Modify”. The sample number/name can be changed manually or the number can be input by a scanner.

Step 6: Remove the test device(test cassette) from the Wellion® INOVO Test Kit and from the foil pouch, and refer to the test kit package insert for the specific test instructions (e.g. sample addition).

Step 7: After preparation of one test cassette, click “Start” on the analyzer display, to start the incubation countdown of the respective test cassette. The analyzer will count down the respective incubation time according to the item (“Project”) you test. Then you can start to sample, prepare and incubate the next test device in the same way. You can complete all samples’ preparation and incubation like this.

Note: After you have clicked the “Start” button for the incubation of the current test cassette, you cannot immediately start the incubation countdown of the next sample on the analyzer (the preparation of the next sample is taken into account as a second and parallel countdown).

Step 8: Insert the test device into the test device slot* sufficient time before the respective incubation countdown is completed. After completion of one countdown, the analyzer will measure the test cassette immediately and, then, display the test result.

***Note:** The test cassette front has to be completely moved inside the analyzer slot, without gap at the front. Carefully check the right direction to move inside which is indicated by an arrow imprint on the test cassettes upper side.


Step 9: Repeat step 6 to 8 to complete testing of all samples.

Insert this direction into the analyzer' cassette slot.

Result

Click the “Result” icon of the main menu interface to enter the following Record interface (Figure 19). You can check the test results and classified statistical operation.

(Figure 19)

Results

The analyzer will automatically save the test result after each completed test; and the user can click the “Result” icon of main menu to check.

1. Select the date range and click the “Confirm” button to view the respective desired records.
2. Use the “Filter” field to select the barcode of tested items, and click “Confirm” to filter out the relevant item records accordingly.
3. Click “Upload” button to upload. Select “All” to upload all records, or select the respective results and click “Select” to upload selected records.
4. Click “Print” button to print results. Select “All” to print all results or select the respective results and click “Select” to print selected records.
5. Click “Export” button to export the selected records or all records.
6. Click “Delete” to delete the selected record or all records.
7. Click “Other” to check detail information of the selected record.

Statistics

Click “Statistics” button to enter the following interface (Figure 20).

(Figure 20)

Select the date range and click the “Confirm” button to view the statistical data of different user account’ results.

Click “Export” to export the statistics results.

Projects

Click the “Projects” icon of the main menu interface to enter the following interface (Figure 21).

(Figure 21)

Ratio

Click the “Ratio” button to enter the following interface (Figure 22).


(Figure 22)

This function should be used by technicians or the ordinary user's who follow the manufacturer's instructions to modify the ITEM's RATIO.

Settings

Click the “Settings” button of the “Projects” interface to enter the following display (Figure 23).

(Figure 23)

1. Users can check the test Item list stored in the analyzer.
Remove item: Select the respective item, click the "Delete" button and click "Confirm" to delete the selected item.
Set item to current: Select the respective item, click "Set" and, then, click "Confirm" to set the test item for the future tests.
2. The function of "Reference" should be used by technicians or the ordinary users who follow the manufacturer's instructions to modify the reference range. Technicians can add new test items or edit the reference range in this interface. The ordinary user need to follow the manufacturer's instructions to add new test items or edit the reference range.

Maintenance and Care

The Wellion® INOVO FIA Analyzer does not require special maintenance or cleaning. Basic maintenance is to keep the analyzer clean.

Maintenance

- Please check the analyzer and accessories regularly.
- Please confirm whether the power outlet is properly grounded. If not, please replace the power outlet.
- Please check whether the power cord is deformed or broken by visual inspection. If there is something wrong with the power cord, please replace by a new one immediately.
- The analyzer uses a built-in thermal printer, which does not require a cartridge change, just require a replacement of the thermal printing paper when needed. When the thermal printing paper has run out, open the built-in thermal printer's cover of the analyzer, and replace the thermal printing paper with a new one, then close the cover.

Protection

- Before cleaning the analyzer, please turn off the power switch and disconnect the power cord.
- Take care to avoid getting liquids, dirt, blood or buffer into the analyzer through the slot or data port.
- Please wipe off the dirt with a soft wet cloth and wipe off the dust with a soft dry cloth when cleaning the housing of the analyzer.

Notes:

- Please do not use gasoline, diluents or other organic solvents to wipe the analyzer; otherwise it will cause the paint shell peeling, damage or deformation.

Troubleshooting

Fault phenomenon	Reason analysis	Solution
Analyzer cannot start	Power switch is not turned on	Turn on the switch
	Power adapter is not connected	Reconnect the power adapter
Screen cannot start	Screen breakdown	Please contact your local dealer or manufacturer
	Operating system problem	
Software system failure	Operating system running fault	Please record the complete error message and contact your local dealer or manufacturer
	Test analysis software cannot start	
	Other error messages appear during software running	
Abnormal sound during testing	The test device holder may be stuck	Turn off the analyzer and turn it on again
	Mechanical motion failure	Please contact your local dealer or manufacturer
Sudden stop during testing	Power interruption	Restart the analyzer and retest
	Communication failure	
	Problem still exists	Please contact your local dealer or manufacturer
Abnormal test result	Abnormal test result	Please contact your local dealer or manufacturer
	Contamination problem	Reduce contamination
Other fault	If another fault occurs	Please contact your local dealer or manufacturer

C line error	Blank test cassette measured	Add sample and buffer
	The test device is inserted in the wrong direction	Insert the test device in the correct direction
	The sample solution has not reached (flowed through) the C line	Wait for the required test incubation time

CHAPTER 6: TECHNICAL INFORMATION

System Specifications

Feature	Specification
Principle	Fluorescence Immunochromatography
Model	Wellion® INOVO FIA Analyzer
Dimensions	215 mm x 303 mm x 159 mm (W x D x H)
Weight	2.5 kg
Operation system	Linux
Channel	Single
Display size	7 inches 24 bits color touch screen
Sample type	Whole blood / Serum / Plasma / Urine / Hair / Swab
Operating temperature	10 °C to 30 °C
Storage temperature	-20 °C to 55 °C
Relative humidity	0 to 93 %
External port	USB / RS-232 / LAN
Connectivity	USB port*, Ethernet port, double serial port (serial port 1: automatic LIS uploading; serial port 2: debugging) *USB port, Use of computer mouse: standard computer mouse can be connected and used to operate analyzer
Power supply	Adapter output: DC 24V 2.5A; Adapter input: 100 - 240 V AC; 50/60 Hz
Printer	Built-in thermal printer

Limitations

Do not place the analyzer in a position that is difficult to operate or inconvenient to disconnect.

Please use the analyzer strictly according to the User's Manual provided by manufacturer.

Do not insert anything into the test device slot except the test device provided by or purchased from the manufacturer.

If the tested samples are potentially contagious, please take protective measures such as wearing protective gloves to avoid the skin contacting the sample.


Please dispose of the analyzer that is no longer used according to the procedures and rules of destruction for Class B electric devices.

Please dispose the used test device according to “Medical Waste Management Regulations” to avoid biohazards.

Dispose the tested samples and materials carefully. Treat them as if they are infectious materials. Follow proper precautions and obey all local regulations when disposing them.

INDEX OF SYMBOLS

These symbols may appear on the packaging and in the instructions of the Wellion® INOVO FIA Analyzer.

Symbols	Explanation
	Consult instructions for use
	Manufacturer
	Date of manufacture
	Temperature limit
	Contains sufficient for <n> tests
	Use-by date
LOT	Batch Code
REF	Catalogue number
#	Model number
SN	Serial number
	Do not re-use
	Caution
	Direct current
	Alternating current
IVD	In vitro diagnostics
	Dispose items according to local relevant laws regarding disposal and recycling
	Keep away from sunlight
	Fragile, handle with care

	Keep dry
	Humidity limit
	Do not use if package is damaged

MANUFACTURER'S WARRANTY

The producer warrants the Wellion® INOVO FIA Analyzer against defects in materials and workmanship for a period of two years from the date of purchase. The warranty is lost if the instrument is misused, poorly maintained or is tampered with.

Liability under this warranty is limited to the repair of defective parts or - at the discretion of the producer - to the replacement of the instrument. The right to rescind the purchase agreement exists only if the replacement is also found to be defective. Claims other than these will not be dealt with.

The warranty is not valid if the damage results from misuse, maltreatment, tampering, human error and the use of extreme.

The warranty period is not to be extended by any claim made under this warranty.

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
The screenshot displays the 'Wellion INOVO FIA Analyzer' software interface. At the top, there is a 'Test' button. Below it, the 'Sample No.' field is empty, while the 'Serial' field contains '211102001'. The 'Sample type' is set to 'Whole blood'. There are buttons for 'Use scanner', 'Info edit', and 'Current user' (admin). The 'Test type' is set to 'Instant test'. Below this, there are fields for 'Sample No.', 'Serial', 'Test project', and 'Sample type'. A 'Read ID' button is positioned to the right. A table with columns 'Result', 'Result', and 'Unit' is visible, with a 'Test' button to its right. A 'Print' button is located at the bottom right. The footer shows the date '11-02-2021', a 'Menu' button, and the time '16:54:46'.

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© MED TRUST, Wellion[®] and INOVO are registered trademarks.