DE Wellion MICRO Teststreifen

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REFERENZEN

1. ADA Clinical Practice Recommendations. 2010
2. According to ISO 15197:2013

EN Wellion MICRO Test strips

PRINCIPE AND INTENDED USE

The Wellion MICRO Blood Glucose Test Strips are thin strips with a chemical reagent system. They work with the PDA of the Wellion MICRO PDA to measure the glucose concentration in whole blood. Blood is applied to the end tip of the test strip. The blood is then automatically absorbed into the reaction cell. The reaction takes place in the reaction cell. A transient electrical current is formed during the reaction which is detected by the meter. The blood glucose concentration is then calculated based on the electrical current. The results are then shown on the meter display. The meters are calibrated to display plasma equivalent results.

For in vitro diagnostic use. Test strips are to be used only outside the body for testing purposes. For self-testing and professional use.

COMPOSITION

Each test strip contains the following reactive chemicals: Glucose Oxidase < 25 IU Mediator < 300 µg

Each test strip vial contains a drying agent.

STORAGE AND HANDLING

[1] Test strips should be stored in their protective vial. The vials cap must be tightly closed. This is to keep the test strips in good working condition [2] Store test strips in a cool, dry place at room temperature, 5-30°C (41-86°F). Store them away from heat and direct sunlight [3] Do not freeze or refrigerate [4] To ensure accurate results, use the test strips at room temperature [5] Keep the colored side of the Wellion MICRO teststrip up and blank white side down when you insert the strip contact bars into the strip port. [6] Do not store the test strips outside of their protective vial. Test strips must be stored in their original vial. The cap must be tightly closed. [7] Do not store or use the test strips in a humid place such as a bathroom [8] Do not store the meter, the test strips or control solution in bleach or cleaners that contain bleach. [9] Do not transfer the test strips to a different tip or vial [10] Place the test strip vial immediately after removing a test strip. [11] Use the test strip immediately after removing it from the vial [12] Do not use your test strips past the unopened expiration date. The expiration date is printed on the vial. Using test strips past the expiration date may produce incorrect test results. Note: All expiration dates are printed in Year-Month format. 2021-01 means end of January 2021

[13] A new vial of test strips may be used for 6 months after first being opened. The opened vial expiration date is 6 months after the date the vial was first opened. Write the opened vial expiration date on the vial label after opening.

PRECAUTIONS

[1] For in vitro diagnostic use, the test strips are to be used only outside the body. The test strips are to be used for testing purposes. [2] Do not use test strips after the expiration date that is shown on the vial. Used expired test strips may give incorrect blood glucose readings. [3] Do not use test strips that are torn, bent, or damaged in any way. Do not freeze or refrigerate. [4] The same method must be applied to the tip of the test strip. Do not apply blood or control solution to the top of the test strip as this may result in an inaccurate reading. [5] Discard the vial and any unused test strips 6 months after you first open it. Constant exposure to air may destroy chemicals in the test strip. This damage can cause incorrect readings. [6] Keep the test strip vial away from children and animals. [7] Consult your healthcare professional before making any changes to your treatment plan.

MATERIALS PROVIDED

- Test Strips
- Package Insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Meter
- Sterile Lancets
- Lancing Device
- Control solution

INSTRUCTIONS FOR USE

See your User's Manual for complete instructions for blood sample collection before use.

1. Open the cap of the test strip vial. Remove a test strip for testing. Replace the cap immediately. This is to protect the remaining test strips from moisture in the air.

2. Run the blood glucose test following the User's Manual.

3. The blood glucose test result will be shown on the meter display window. This result should fall within the target range. Your healthcare professional should recommend your target range. If your blood glucose test results are higher or lower, always consult your healthcare professional before making any changes to your treatment plan.

IMPORTANT: Wellion MICRO Blood Glucose Monitoring Systems allow alternative site testing for forearm and palm testing in addition to fingertip testing. There are important differences between forearm palm and fingertip samples that you should know. Important information about forearm and palm glucose testing:

[1] When blood levels are changing rapidly such as after a meal, insulin dose or exercise, blood from the fingertips may show these changes more rapidly than blood from other areas.

[2] Fingertips should be used if testing is within 2 hours of a meal, insulin dose or exercise and any time you feel glucose levels are changing rapidly.

[3] You should test with the fingertips anytime there is a concern for hypoglycemia or you suffer from hypoglycemia unawareness.

RANGE OF EXPECTED VALUES

Blood glucose monitoring requires the help of a healthcare professional. Together you can set your own range of expected blood glucose values, arrange your testing times and discuss the meaning of your blood glucose results. Expected blood glucose levels for people without diabetes:

Time Range, mg/dL Range, mmol/L

Fasting and before meals 70-100 3,9-5,6

2 Hours After Meal Less than 140 Less than 7,8

CHECKING THE SYSTEM

Your blood glucose meter must be handled carefully. See your User's Manual for detailed instructions for meter care. The quality control test should be used to check that the meter and test strips are working together properly. Follow the test procedure in your User's Manual to run a quality control test. Two ranges, CTRL 1 and CTRL 2 are shown on the test strip vial label. Control Solution 1 is sufficient for most all self-testing needs. If you think your meter or strips are not working correctly, you may also want to do a test with your dealer for help.

For confirmation of results, Control Solution 1 tests should fall within the CTRL 1 range, and Control Solution 2 test should fall within the CTRL 2 range. When testing with Control Solution 1, make sure you are matching the results to the CTRL 1 range on the vial label.

CAUTION: If your quality control test results fall outside the control range shown on the test strip vial, DO NOT use the system to test your blood, as the system may not be working properly. If you cannot correct the problem, contact your dealer for help.

LIMITATIONS

[1] The Wellion MICRO test strips and the Wellion MICRO PDA have been designed, tested and proven to work together effectively to provide accurate blood glucose measurements. Do not use components from other brands.

[2] The Wellion MICRO Test Strips are for testing with fresh capillary whole blood. Do not use with serum or plasma samples. [3] Do not use for testing newborns. [4] Very high levels of triglycerides (triglycerides up to 950 mg/dL, 13,0 mmol/L) or cholesterol up to 13 mmol/L (150 mg/dL) may have no major effect on blood glucose test results.

[5] The Wellion MICRO Blood Glucose Monitoring Systems show to work properly in studies at altitudes up to 3,048 meters. [6] Severely ill persons should not run the glucose test with the Wellion MICRO Blood Glucose Monitoring Systems [10] Blood samples from patients in shock, or with severe dehydration or from patients in a hyperosmolar state (with or without ketosis) have not been tested and are not recommended for testing with Wellion MICRO Blood Glucose Monitoring System. [11] Disposal of blood samples as they are infected materials. Follow proper precautions when disposing of materials.

PERFORMANCE CHARACTERISTICS

Reproducibility, Precision